

## Autofil™ Centrifuge Bottles

RLS Part No.	Description	Dimensions		Max. RCF (g)*
		Inches	mm	
1801-RLS	250ml Polypropylene bottle with Scew cap, nonsterile, 36 pk	2.32x5.01	59x129	8250
1803-RLS	250ml Polypropylene bottle with Sealing cap, nonsterile, 36 pk	2.32x5.01	59x129	27500
1805-RLS	500ml Polypropylene bottle with Scew cap, nonsterile, 24 pk	2.68x6.61	68x168	8250
1807-RLS	500ml Polypropylene bottle with Sealing cap, nonsterile, 24 pk	2.68x6.61	68x168	13500
1802-RLS	250ml Polycarbonate bottle with Scew cap, nonsterile, 36 pk	2.32x5.01	59x129	8250
1804-RLS	250ml Polycarbonate bottle with Sealing cap, nonsterile, 36 pk	2.32x5.01	59x129	31700
1806-RLS	500ml Polycarbonate bottle with Scew cap, nonsterile, 24 pk	2.68x6.61	68x168	8250
1808-RLS	500ml Polycarbonate bottle with Sealing cap, nonsterile, 24 pk	2.68x6.61	68x168	13850

\* bottles must be filled 80% of capacity.

RLS Part No.	Description	Suggested List Price	Dealer Price		
			(40-99 cases) 30% Discount	(100-199 cases) 40% Discount	(200+ cases) 50% Discount
1801-RLS	250ml Polypropylene bottle with Scew cap, nonsterile, 36 pk	\$258.55	\$180.99	\$155.13	\$129.28
1803-RLS	250ml Polypropylene bottle with Sealing cap, nonsterile, 36 pk	\$639.25	\$447.48	\$383.55	\$319.63
1805-RLS	500ml Polypropylene bottle with Scew cap, nonsterile, 24 pk	\$242.25	\$169.58	\$145.35	\$121.13
1807-RLS	500ml Polypropylene bottle with Sealing cap, nonsterile, 24 pk	\$560.76	\$392.53	\$336.46	\$280.38
1802-RLS	250ml Polycarbonate bottle with Scew cap, nonsterile, 36 pk	\$435.25	\$304.68	\$261.15	\$217.63
1804-RLS	250ml Polycarbonate bottle with Sealing cap, nonsterile, 36 pk	\$732.49	\$512.74	\$439.49	\$366.25
1806-RLS	500ml Polycarbonate bottle with Scew cap, nonsterile, 24 pk	\$460.50	\$322.35	\$276.30	\$230.25
1808-RLS	500ml Polycarbonate bottle with Sealing cap, nonsterile, 24 pk	\$598.99	\$419.29	\$359.39	\$299.50

\* Relative Centrifugal Force